


Product Datasheet

SYP antibody (orb66899)

Description	Rabbit polyclonal to SYP. Synaptophysin (SYP) is an integral membrane protein of small synaptic vesic...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SYP
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	34
Uniprot ID	P08247
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...