



Product Datasheet

Enkephalin antibody (orb6690)


Description	Enkephalin antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Enkephalin. Please contact us for the exact immunogen sequence. The peptide is available as orb12923.
Target	Enkephalin
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	31 kDa
Uniprot ID	P22005, P04094
NCBI	001129162
Dilution Range	WB: 1:100-1500, IHC-P: 1:100-500, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



IHC-P image of rat testis tissue using a...



Western blot analysis of mouse brain (la...



Immunofluorescence analysis of rat testi...