

Product Datasheet

Human LILRB5/CD85c/LIR-8 ELISA Kit (orb669045)

Description	Human LILRB5/CD85c/LIR-8 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: G24-G458
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	O75023
Sensitivity	15 pg/ml
Expiration Date	Please enquire.

