


Product Datasheet

AKAP12 Antibody (orb669056)


Description

AKAP12 Antibody


Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived mouse AKAP12 recombinant protein (Position: E387-H1620).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	300 kDa
Uniprot ID	Q9WTQ5
Dilution Range	Western blot, 0.1-0.25µg/ml, Mouse Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Mouse Immunocytochemistry/Immunofluorescence, 5µg/ml, Mouse Flow Cytometry, 1-3µg/1x10 ⁶ cells, Mouse, Rat Direct ELISA, 0.1-0.5µg/ml, Mouse
Expiration Date	12 months from date of receipt.



Western blot analysis of AKAP12 using an...



IHC analysis of AKAP12 using anti-AKAP12...



Flow Cytometry analysis of HEPA1-6 cells...