


Product Datasheet

KDM5A antibody (orb66908)

Description	Rabbit polyclonal to KDM5A. Lysine (K)-specific demethylase 5A (KDM5A) is a ubiquitously expressed nu...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human KDM5A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	192
Uniprot ID	P29375
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...