


Product Datasheet

TDT antibody (orb66909)

Description	Rabbit polyclonal to TDT. TDT is a template-independent DNA polymerase that belongs to the DNA polyme...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human TDT
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	58
Uniprot ID	P04053
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...