


Product Datasheet

UBE2B antibody (orb66911)

Description	Rabbit polyclonal to UBE2B. The modification of proteins with ubiquitin is an important cellular mech...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human UBE2B
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	17
Uniprot ID	P63146
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...