

Product Datasheet

MLL antibody (orb66913)

Description	Rabbit polyclonal antibody to MLL
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MLL
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	432
Uniprot ID	Q03164
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.

