


Product Datasheet

BLM antibody (orb66916)

Description	Rabbit polyclonal to BLM. BLM (Bloom syndrome protein) is a magnesium-dependent and ATP-dependent DNA...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human BLM
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	159
Uniprot ID	P54132
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...