

## Product Datasheet

### CDH15 antibody (orb66917)

|                            |                                     |
|----------------------------|-------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to CDH15 |
| <b>Species/Host</b>        | Rabbit                              |
| <b>Reactivity</b>          | Human                               |
| <b>Conjugation</b>         | Unconjugated                        |
| <b>Tested Applications</b> | WB                                  |
| <b>Immunogen</b>           | N-terminal region of human CDH15    |
| <b>Preservatives</b>       | 0.01% NaN <sub>3</sub>              |
| <b>Storage</b>             | -20°C                               |
| <b>Note</b>                | For research use only               |
| <b>Clonality</b>           | Polyclonal                          |
| <b>Purity</b>              | Affinity purification               |
| <b>MW</b>                  | 89                                  |
| <b>Uniprot ID</b>          | <b>P55291</b>                       |
| <b>Dilution Range</b>      | WB:1:1000-1:5000                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.     |

