

Product Datasheet

HOXD11 Antibody (orb669180)

Description
HOXD11 Antibody
Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, ICC, IF, IHC, WB

Applications
Immunogen

A synthetic peptide corresponding to a sequence of human HOXD11 (VREVAFRDYGLERAKWPYR).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


45 kDa

Uniprot ID
P31277
Dilution Range


Western blot, 0.25-0.5µg/ml, Human, Mouse
Immunohistochemistry (Paraffin-embedded Section), 2-5µg/m,
Human Immunocytochemistry/Immunofluorescence, 5µg/ml,
Human Flow Cytometry, 1-3µg/1x10⁶ cells, Human, Mouse,
Rat

Expiration Date


12 months from date of receipt.



Western blot analysis of HOXD11 using anti-HOXD11...



IHC analysis of HOXD11 using anti-HOXD11...



IHC analysis of HOXD11 using anti-HOXD11...