

Product Datasheet

ZHX2 antibody (orb66919)

Description	Rabbit polyclonal antibody to ZHX2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human ZHX2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	92
Uniprot ID	Q9Y6X8
Dilution Range	WB:1:500-1:1000
Expiration Date	12 months from date of receipt.

