

Product Datasheet

YTHDF2 Antibody (orb669214)

Description

YTHDF2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, ICC, IF, IHC, WB

Applications
Immunogen

A synthetic peptide corresponding to a sequence of human YTHDF2 (NDDDFEYLSPQARPNNAYTAMSD).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


65 kDa

Uniprot ID
Q9Y5A9
Dilution Range


Western blot, 0.1-0.25µg/ml, Human, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Human, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 2µg/ml, Human
Flow Cytometry, 1-3µg/1x10⁶ cells, Human, Mouse, Rat

Expiration Date

12 months from date of receipt.



Western blot analysis of YTHDF2 using anti-YTHDF2...



IHC analysis of YTHDF2 using anti-YTHDF2...



IHC analysis of YTHDF2 using anti-YTHDF2...