

Product Datasheet

SARS-CoV-2 (COVID-19) NSP5A Protein, His Tag (orb669262)

Description	0
Reactivity	Mouse
Form/Appearance	Lyophilized
Storage	The product is shipped at ambient temperature. Upon receipt, store it immediately at -20°C for 6 months under sterile conditions. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Protein Sequence	6xHis tag at C-terminal Accession #: YP_009725301 SGFRKMAFPSGKVEGCMVQVTCGTTTTLNGWLDDVVYCPRHVICTSEDMLNPN YEDLLIRKSNHNFLVQAGNVQLRVIGHSMQNCVLKLVDTANPKTPKYKFVRIQP GQTFSVLACYNGSPSGVYQCAMRPNFTIKGSFLNGSCGSGVGFNIDYDCVSFCYM HHMELPTGVHAGTDLEGNFYGPVDRQTAQAAGTDTTITVNVLAWLAAVING DRWFLNRFTTTTLNDFNLVAMKYNIEPLTQDHVDILGPLSAQTGIAVLDMCASLK ELLQNGMNGRTILGSALLEDEFTPFDVVRQCSGVTFQ
Source	E. coli
MW	34.6KD
Expiration Date	6 months from date of receipt.

