


Product Datasheet

SMARCC1 antibody (orb66928)

Description	Rabbit polyclonal to SMARCC1.SWI/SNF complex subunit SMARCC1 is a member of the SWI/SNF family prote...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SMARCC1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	123
Uniprot ID	Q92922
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...