

Product Datasheet

Biotinylated Human TNFR2 / CD120b / TNFRSF1B Protein (orb669285)

Description

Biotinylated Human TNFR2, His, is

(HEK293). It contains AA Leu 23 - A...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

TNFR2

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His & C-Avi

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 28.7 kDa. The protein migrates as 40-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

AAA36755

Purity

90%

Source

Biotinylated Human TNFR2, His, (orb669285) is expressed from human 293 cells (HEK293). It contains AA Leu 23 - Asp 257 (Accession # P20333-1).

MW

28.7 kDa

Uniprot ID

P20333

NCBI

AAA36755

Expiration Date

6 months from date of receipt.