

Product Datasheet

Human ILDR2 Protein (orb669286)

Description	Human ILDR2, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Glu 186 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	ILDR2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 45.1 kDa. The protein migrates as 45-50 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.
Protein Sequence	NP_955383.1
Purity	95%
Source	Human ILDR2, Fc Tag (orb669286) is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Glu 186 (Accession # Q71H61-1).
MW	45.1 kDa
Uniprot ID	Q71H61
NCBI	NP_955383.1
Expiration Date	6 months from date of receipt.