

Product Datasheet

**SARS-CoV-2 (COVID-19) S2 protein [Out of stock]
(orb669291)**

Description

SARS-CoV-2 S2 protein, His Tag

human 293 cells (HEK293).

Reactivity Virus

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Spike S2

Preservatives 20 mM PB, 300 mM NaCl, pH7.4

Form/Appearance Liquid

Storage -70°C

Tag C-10×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 60.0 kDa. The protein migrates as 65-115 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence QHD43416.1

Purity ● The product MUST be stored at -70°C or lower upon receipt; ● -70°C for 3 months under sterile conditions.

Source SARS-CoV-2 S2 protein, His Tag (orb669291) is expressed from human 293 cells (HEK293).

MW 60.0 kDa

NCBI [QHD43416.1](#)

Expiration Date 6 months from date of receipt.