

Product Datasheet

SARS-CoV-2 (COVID-19) S protein RBD [Out of stock] (orb669298)

Description	It's been reported that Coronavirus can infect the human respiratory epithelial cells through intera...
Reactivity	Virus
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Spike RBD
Preservatives	10 mM PB, pH7.4
Form/Appearance	Liquid
Storage	-70°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 26.5 kDa. The protein migrates as 32 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	QHD43416.1
Purity	95%
Source	SARS-CoV-2 S protein RBD, His Tag (orb669298) is expressed from human 293 cells (HEK293).
MW	26.5 kDa
NCBI	QHD43416.1
Expiration Date	6 months from date of receipt.