

Product Datasheet

Rhesus macaque Properdin Protein (orb669299)

Description

Properdin is also known as Complement

the alternate path...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Properdin

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 50.6 kDa. The protein migrates as 55-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Rhesus macaque Properdin, His Tag (orb669299) is expressed from human 293 cells (HEK293). It contains AA Asp 28 - Leu 469 (Accession # F7CGG5-1).

MW

50.6 kDa

NCBI

[XP_014982784.1](#)

Expiration Date

6 months from date of receipt.