

Product Datasheet

SARS-CoV-2 Nucleocapsid (NP) protein (orb669317)

Description	Recombinant SARS-CoV-2 Nucleocapsid Protein is produced by E
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 um filtered solution of PB, pH 7.4
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PB, pH 7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at - 20°C for 3 months.
Tag	N-6His
Note	For research use only
Protein Sequence	1-419
Source	E. coli
Uniprot ID	P0DTC9
Expiration Date	6 months from date of receipt.