

Product Datasheet

SARS-CoV-2 Spike (RBD) protein (orb669321)

| | |
|-------------------------|---|
| Description | Recombinant SARS-CoV-2 Spike Protein is produced by Insect |
| Endotoxins | Less than 0.1 ng/ug (1 EU/ug) |
| Preservatives | Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4 |
| Form/Appearance | Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4 |
| Storage | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| Tag | C-8His |
| Note | For research use only |
| Protein Sequence | 319-532 |
| Source | Insect Cells |
| Uniprot ID | P0DTC2 |
| Expiration Date | 6 months from date of receipt. |