

Product Datasheet

SARS-CoV-2 Spike (S1+S2 ECD) protein (orb669323)

Description	Recombinant SARS-CoV-2 Spike Protein is produced by Insect
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 um filtered solution of PBS, pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze- thaw cycles.
Tag	C-8His
Note	For research use only
Protein Sequence	15-1214
Source	Insect Cells
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.