

Product Datasheet

Human Acyl-Protein Thioesterase 1 protein (orb669414)

Description	Recombinant Human Acyl-Protein Thioesterase 1 is produced by our E
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Protein Sequence	Met1-Asp230
Source	E.coli
Uniprot ID	O75608
Expiration Date	6 months from date of receipt.