

Product Datasheet

Human Fructose-Bisphosphate Aldolase A protein (orb669428)

Description	Recombinant Human Fructose-Bisphosphate Aldolase A is produced by our E
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 20% Glycerol, pH 8.0
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, 20% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Pro2-Tyr364
Source	E.coli
Uniprot ID	P04075
Expiration Date	6 months from date of receipt.