

Product Datasheet

Human Phosphopantothenoylcysteine Decarboxylase protein (orb669435)

Description	Recombinant Human PPC-DC is produced by our E
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 50mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 50mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Protein Sequence	Met1-Ser204
Source	E.coli
Uniprot ID	Q96CD2
Expiration Date	6 months from date of receipt.