

Product Datasheet

CNN1 antibody (orb66945)

Description	Rabbit polyclonal to CNN1. Calponin-1 (CNN1) is a thin filament-associated protein implicated in the ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human CNN1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	33
Uniprot ID	P51911
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.

