

Product Datasheet

Human XPNPEP2 protein (orb669468)

Description	Recombinant Human Aminopeptidase P2 is produced by our Mammalian expression system and the target ge...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	His22-Ala650
Source	Human Cells
Uniprot ID	O43895
Expiration Date	6 months from date of receipt.