

Product Datasheet

Human MINPP1 protein (orb669471)

Description	Recombinant Human MINPP1 is produced by our Mammalian expression system and the target gene encoding...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Ser31-Leu487
Source	Human Cells
Uniprot ID	Q9UNW1
Expiration Date	6 months from date of receipt.