

Product Datasheet

Human ALPL protein (orb669479)

Description	Recombinant Human Alkaline Phosphatase, Tissue-Nonspecific Isozyme is produced by our Mammalianexpre...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM HEPES, 150mM NaCl, 2mM MgSO4, 0.1mM ZnCl2, pH 7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM HEPES, 150mM NaCl, 2mM MgSO4, 0.1mM ZnCl2, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Leu18-Ser502
Source	Human Cells
Uniprot ID	P05186
Expiration Date	6 months from date of receipt.