


Product Datasheet

Adiponectin antibody (orb66949)

Description	Rabbit polyclonal antibody to Adiponectin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human ADIPOQ
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	26
Uniprot ID	P52566
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...