

Product Datasheet

Human Glucose-6-Phosphate 1-Dehydrogenase protein (orb669502)

Description	Recombinant Human G6PD is produced by our Mammalian expression system and the target gene encoding A...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Ala2-Leu515
Source	Human Cells
Uniprot ID	P11413
Expiration Date	6 months from date of receipt.