


Product Datasheet

ABL1 antibody (orb66952)

Description	Rabbit polyclonal to ABL1.The ABL1 protooncogene encodes a cytoplasmic and nuclear protein tyrosine ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human ABL1
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	123
Uniprot ID	P00519
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...