


Product Datasheet

CUL4A antibody (orb66954)

Description	Rabbit polyclonal to CUL4A. CUL4A is the ubiquitin ligase component of a multimeric complex involved ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human CUL4A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	88
Uniprot ID	Q13619
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...