

Product Datasheet

Human ADP-Ribosylarginine Hydrolase protein (orb669543)

Description	Recombinant Human ADP-Ribosylarginine Hydrolase is produced by our E
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mMtris,50mMNaCl,1mMDTT,20%Glycerol,PH8.0
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mMtris,50mMNaCl,1mMDTT,20%Glycerol,PH8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Protein Sequence	Met1-Leu357
Source	E.coli
Uniprot ID	P54922
Expiration Date	6 months from date of receipt.