

Product Datasheet

Human Cytoplasmic NCK1 protein (orb669552)

Description	Recombinant Human NCK1 is produced by our E
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, 5% Mannitol, pH 8.0
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 5% Mannitol, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Protein Sequence	Met1-Ser377
Source	E.coli
Uniprot ID	P16333
Expiration Date	6 months from date of receipt.