

## Product Datasheet

### Human SMAD Family Member 4 protein (orb669553)

<b>Description</b>	Recombinant Human SMAD4 is produced by our E
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, pH 8.0
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, pH 8.0.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Met1-Asp552
<b>Source</b>	E.coli
<b>Uniprot ID</b>	<b>Q13485</b>
<b>Expiration Date</b>	6 months from date of receipt.