

Product Datasheet

Human Ubiquitin-Conjugating Enzyme Variant 2 protein (orb669561)

Description	Recombinant Human Ubiquitin-Conjugating Enzyme E2 Variant 2 is produced by our E
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 50mm HEPES,150mM NaCl, pH 7.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mm HEPES, 150mM NaCl, pH 7.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Protein Sequence	Met1-Asn145
Source	E.coli
Uniprot ID	Q15819
Expiration Date	6 months from date of receipt.