


Product Datasheet

ELL2 antibody (orb66958)

Description	Rabbit polyclonal to ELL2.ELL2 belongs to the ELL/occludin family. ELL2 is an elongation factor that...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ELL2 [distinct from the antigen for R0211-1, ELL2 (C1) antibody]
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	72
Uniprot ID	O00472
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...