


Product Datasheet

ELL2 antibody pair (orb66959)

Description	Rabbit polyclonal to ELL2. 50 μl ELL2 (C-term1) &50 μl ELL2 (C-term2).vPair antibodies represe...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	-
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	72
Uniprot ID	O00472
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...