


Product Datasheet

ATF1 antibody pair (orb66963)

Description	Rabbit polyclonal to ATF1. 50 μl ATF1 (N-term) &50 μl ATF1 (C-term) .vPair antibodies represen...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	-
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	29
Uniprot ID	P18846
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...