

Product Datasheet

SARS-CoV-2 (2019-nCoV) 3C-likease protein (orb669653)

Description	Recombinant SARS-CoV-2 (2019-nCoV) 3C-like Proteinase is produced by our E
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of PBS, 1mM EDTA, 10% Glycerol, pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 1mMEDTA, 10%glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Protein Sequence	Ser1-Gln306
Source	E.coli
Uniprot ID	YP_009725301.1
Expiration Date	6 months from date of receipt.