

Product Datasheet

Human Thrombomodulin protein (orb669672)

Description	Recombinant Human Thrombomodulin is produced by our expression system and the target gene encoding A...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Ala19-Ser515
Source	Human Cells
Uniprot ID	P07204
Expiration Date	6 months from date of receipt.