

Product Datasheet

Biotinylated Human ACE-2 protein (orb669689)

Description	Biotinylated Recombinant Human Angiotensin-Converting Enzyme 2 is produced by our Mammalian expressi...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Conjugation	Biotin
Preservatives	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 300mM NaCl, 1mM ZnCl2, 10% Glycerol, pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mMTris-HCl,300mMNaCl,1mMZnCl2,10%Glycerol,PH7.4 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-AVI&His
Note	For research use only
Protein Sequence	Gln18-Ser740
Source	Human Cells
Uniprot ID	Q9BYF1
Expiration Date	6 months from date of receipt.