

Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid protein (orb669691)

Description	Recombinant SARS-CoV-2 (2019-nCoV) Nucleocapsid Protein is produced by our E
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µM filtered solution of 20 mM Tris-HCl, 300 mM NaCl, 1mM EDTA, pH 8.0.
Form/Appearance	Supplied as a 0.2 µM filtered solution of 20 mM Tris-HCl, 300 mM NaCl, 1mM EDTA, pH 8.0.
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Protein Sequence	Met1-Ala419
Source	E. coli
Uniprot ID	QHD43423.2
Expiration Date	6 months from date of receipt.