

Product Datasheet

SARS-CoV-2 (2019-nCoV) Envelope protein (orb669692)

Description	Recombinant SARS-CoV-2 (2019-nCoV) Envelope Protein is produced by our E
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µM filtered solution of 20mM Tris-HCl, 200mM NaCl, pH 8.0.
Form/Appearance	Supplied as a 0.2 µM filtered solution of 20mM Tris-HCl, 200mM NaCl, pH 8.0.
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	N-BBP&His
Note	For research use only
Protein Sequence	Met1-Val75
Source	E. coli
Uniprot ID	QHD43418.1
Expiration Date	6 months from date of receipt.