

Product Datasheet

SARS-CoV-2 (2019-nCoV) S1 protein (orb669694)

Description	Recombinant SARS-CoV-2 (2019-nCoV) S1 Protein is produced by our Mammalian expression system and the...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µM filtered solution of PBS, pH 7.4.
Form/Appearance	Supplied as a 0.2 µM filtered solution of PBS, pH 7.4.
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-Fc
Note	For research use only
Protein Sequence	Gln14-Arg685
Source	Human Cells
Uniprot ID	QHD43416.1
Expiration Date	6 months from date of receipt.