

Product Datasheet

2019-nCoV S-trimer protein (orb669704)

Description	Recombinant 2019-nCoV S-trimer Protein is produced by our Mammalian expression system and the target...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µM filtered solution of PBS, pH 7.4.
Form/Appearance	Supplied as a 0.2 µM filtered solution of PBS, pH 7.4.
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Cys15-Gln1208(D614G)
Source	Human Cells
Uniprot ID	QHD43416.1
Expiration Date	6 months from date of receipt.