

Product Datasheet

SARS-CoV-2 S-Protein ACE2 Binding Domain Antibody (orb669721)

Description

Our SARS-CoV-2 S-Protein ACE2

	from Biorbyt is mouse monoclonal. It d...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant SARS-CoV-2 S-Protein construct containing the ACE2 binding domain expressed in and purified from E. coli.
Target	SARS-CoV-2 S-Protein ACE2 Binding Domain
Preservatives	PBS + 50% glycerol and 5 mm NaN ₃
Form/Appearance	Protein G Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	25 or 142 full-length
Uniprot ID	P0DTC2
Dilution Range	WB: 1:2000, ICC: 1:1000
Expiration Date	12 months from date of receipt.