


Product Datasheet

MYLIP antibody (orb66979)

Description	Rabbit polyclonal to MYLIP. Myosin regulatory light chain interacting protein (MYLIP) is also known a...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MYLIP
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	50
Uniprot ID	Q8WY64
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...